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- FDA Approved
- US Patented
- International Patented (CA, EA, Asia, etc.)
- Short Pulse RF

- No Side Effect
- Less Downtime
- Less Bleeding
- Safe for All Kinds of Skin Type

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


THE DIFFERENCE BETWEEN NON INSULATED ELECTRODES AND INSULATED ELECTRODES

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

PMCID : PMC4643287

Published online 2015 Nov 13, doi:10.1038/srclp16735

Electromagnetic Initiation and Propagation of Bipolar Radiofrequency Tissue Reactions via Invasive Non-Insulated Microneedle Electrodes

Jongju Na, Zhenlong Zheng, Christopher Dannaker, Sang Eun Lee, Jin-Soo Kang and Sung Bin Cho Author information / Artuclue bites / Copyright and License information

Abstract

Radiofrequency (RF) energy can be emitted into the skin, either non-or invasively, via a monopolar mode that utilizes an active electrode and a grounded electrode or via a bipolar mode that employs two active electrodes. In this experimental study of RF tissue reactions, bipolar RF energy was emitted in vivo to micropig skin at varying microneedle penetration depths, signal amplitudes, and conduction times. Immediately after RF treatment, skin samples exhibited RF-induced coagulation columns of thermal injury, separately generated around each microneedle in the dermis. In ex vivo bovine liver tissue, the thermal coagulation columns were found to be concentrated maximally around the pointed tips of each electrode. After a RF conduction time of 2 seconds, the individual areas of thermal coagulation began to converge with neighboring RF-induced coagulation columns; the convergence of coagulation columns was found to start from the tips of neighboring electrodes.

Radiofrequency (RF) refers to high frequency alternating electrical current at the frequency range traditionally used for radio-wave communication. Electromagnetic signal, including RF, induces an electrothermal reaction in targeted tissues, the patterns of which depend on the resistance of the tissue. RF energy can be emitted into the skin either non-or invasively via a monopolar mode that utilizes an active electrode and a grounded electrode or via a bipolar mode that employs two active electrodes. In monopolar modes, an electrical circuit formed by an electron current that flows from the active electrode to the grounded electrode is generated in the patient's body. Meanwhile, bipolar electrosurgery systems induce an electrical circuit between the two active electrodes that is limited to regionally targeted tissues. In the skin these active electrodes emit an electron current that flows through the shortest path in the target tissue between the electrodes, theoretically limiting the depth of the thermal response induced by the electromagnetic energy. Invasive RF systems, which deliver electromagnetic energy through electrodes that penetrate into target tissues, offer advantages of a deeper layer of treatment in a non-contiguous pattern, compared with noninvasive systems. In a recent study, a monopolar 0.4-MHz RF system equipped with a single penetrating electrode induced thermal coagulation in ex vivo bovine liver tissue that started from the tips of non-insulated electrodes and formed a rim of coagulated tissue around the entire length of the needles with increasing energy levels. Therein, the thickest rim of coagulated tissue was formed around the tips of the electrodes, suggesting that non-insulated penetrating electrodes can be used to effectively and safely deliver RF energy to the skin while preserving the epidermis.

Hyun Jun Park

Published by KSLMS

Park, HJ, Kim HM, Oh MJ. Clinical Study of Facial Wrinkle Treatment with Fractional Microneedle Radio Frequency System / Med LASER 44:631-636(2012)

Kyu Young Seo

Published by Lasers in Surgery and Medicine

Seo KY, Kim DH, Yoon MS Et al. Skin Rejuvenation by Microneedle Fractional Radiofrequency Treatment in Asian Skin; Clinical and Histological Analysis / Med LASER 2014;3(2):59-64

Sang Ju Lee

American Society for Dermatologic Surgery

Lee, SJ, et al. Use of Fractionated Microneedle Radiofrequency for the Treatment of Inflammatory Acne Vulgaris in 18 Korean Patients / Dermatol Surg 2012;1-6

Ji Hun Park

Published by Laser Surg Med / 2013

Park JH, Shin JU, Cho SH Et al. The use of microneedle fractional radiofrequency system in wrinkle reduction and skin tightening

Kyu Young Seo

Published by KSLMS

Na, JJ, Et al. The Efficacy and Safety of Bipolar Radiofrequency Treatment with NonInsulated Penetrating Microneedles for Acne Vulgaris and Acne Scars / Med Laser 2015;4(1):10-15

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©2023 ViOL Co., Ltd. All rights reserved.